



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Advanced Laboratory Range

Advanced Laboratory Freezer RAFR21263

The Labcold Advanced ranges of refrigerators and freezers represent the highest levels of performance, constructional quality, robustness, and appearance available in the market place.

Heavy duty all-metal construction, with optional stainless steel exterior provides strength and durability with ease of cleaning or sanitisation.

Microprocessor control, together with ducted air flow and advanced auto defrost system provide class leading temperature stability and uniformity.

Options available: Sample storage racking; Remote alarm contacts; Access ports; Internal power sockets with safety trip; Over and under temperature safety trips; Temperature recorders/loggers.

Advanced models may also be customised for specific applications – please ask us.



Specifications

Product Code	RAFR21263
Gross Capacity	650 litres
Style	Free standing
UK Warranty	2 Years
External Dimensions (mm) HxWxD	1990mm x 710mm x 800mm
Net Weight	130kg
Shelves	4
Maximum Shelf Loading	70kg

Features

- Microprocessor temperature control
- Audible and visual alarm
- Temperature uniformity
- 2 year parts and labour warranty (UK only)
- Digital temperature display
- Optimised air circulation
- Automatic defrosting

Models

	Model	Gross Capacity	Door Style	Lock Type
	Advanced Laboratory Freezer RAFR21263	650 litres	Solid	Key lock
	Advanced Laboratory Freezer RAFR44263	1300 litres	Solid	2 key locks