



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree ATEX Refrigerators and Freezers

## Sparkfree ATEX Fridge RLPRO217-ATEX

This 66 litre capacity sparkfree ATEX refrigerator is manufactured with an interior which is licensed as a Zone hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

Despite being compact, this refrigerator is fitted with a door lock, digital temperature display with alarms and a condensate collection vessel and backed by a free 2 year UK parts and labour warranty.

WARNING: Do not store open containers of volatile substances in this fridge

### Specifications

Product Code	RLPRO217-ATEX
Gross Capacity	66 litres
Style	Bench top
UK Warranty	Two years
External Dimensions (mm) HxWxD	740 x 450 x 510
Internal Dimensions (mm) HxWxD	590 x 333 x 300
Net Weight	34kg
Shelves	3

## Specifications (cont.)

Temperature Range	0°C to 10°C
-------------------	-------------

Internal Material	Painted Aluminum
-------------------	------------------

Door Style	Solid
------------	-------

Temperature Display	Yes
---------------------	-----

Door Open Alarm	Yes
-----------------	-----

Condensate Collection	Yes
-----------------------	-----

Manual Defrost	No
----------------	----

Eco Friendly Refrigerant	Yes
--------------------------	-----

Factory Set Point	5°C
-------------------	-----

External Material	Powder Coated Steel
-------------------	---------------------

Lock Type	Key lock
-----------	----------

High/Low Temperature Alarm	Yes
----------------------------	-----

Door Earth Connection	Yes
-----------------------	-----

Off Cycle Auto Defrost	Yes
------------------------	-----



Power Consumption Kwh/24hrs	0.6
--------------------------------	-----

Ambient Operating Temperature	16°C to 32°C
----------------------------------	--------------

## Features

- Digital temperature display
- Door open alarm
- Door earth connection
- Off cycle auto defrost
- High/low temperature alarm
- Door lock
- Condensate collection
- 2 years UK parts and labour warranty

## Models

	Model	Gross Capacity	Door Style	Lock Type
	Sparkfree ATEX Fridge RLPRO217-ATEX	66 litres	Solid	Key lock
	Sparkfree ATEX Freezer RLVFO217-ATEX	66 litres	Solid	Key lock