



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree ATEX Refrigerators and Freezers

Sparkfree ATEX Freezer RLVF0417-ATEX

This 124 litre capacity sparkfree ATEX freezer is manufactured with an interior which is licensed as a Zone hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

This under counter laboratory freezer is also fitted with a door lock, digital temperature display with alarms and a condensate collection vessel and backed by a free 2 year UK parts and labour warranty.

WARNING: Do not store open containers of volatile substances in this freezer



Specifications

Product Code	RLVF0417-ATEX
Gross Capacity	124 litres
Style	Under counter
UK Warranty	Two years
External Dimensions (mm) HxWxD	830 x 600 x 600
Internal Dimensions (mm) HxWxD	680 x 460 x 410
Temperature Range	-18°C to -25°C
Factory Set Point	-20°C

Features

- Digital temperature display
- Door open alarm
- Door earth connection
- High/low temperature alarm
- Door lock
- 2 years UK parts and labour warranty

Models

	Model	Gross Capacity	Door Style	Lock Type
	Sparkfree ATEX Freezer RLVF0417-ATEX	124 litres	Solid	Key lock
	Sparkfree ATEX Fridge RLPR0517-ATEX	150 litres	Solid	Key lock

Accessories



Labcold Stacking Kit